

Test schedule

Follow the schedule to monitor your washer for 14 working days. Use alongside your standard load.

WASHER ID:

DAY / DATE	LOADING POSITION	CYCLE NUMBER	PROTEIN QTY [µg]			
1	1, 2, 3, 4					
2	1					
3	2					
4	4					
5	1					
6	2					
7	4					
8	1, 2, 3, 4					
9	1					
10	2					
11	4					
12	1					
13	2					
14	4					

Ordering information

Product Code	Product Name	Packaging
018.1000	VeriTest Tag, Stainless Steel with BB Test Soil	Pack of 20
018.1010	VeriTest Multi Block	Single Unit
018.2000	VeriTest Blue	Pack of 20